



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipe Friction Training Panel

Product Code :
LBNY-0005-12500041



Description :

The trainer is supplied from the water mains. The flow rate is adjusted using valves and measured using a rotameter. Four different pipe sections are fitted to a panel. Each section can be individually selected using ball valves. The measuring points are 1m apart. The pressure drop over a pipe section is measured using a digital, battery-operated differential manometer. To keep measuring inaccuracies as low as possible, the pressure is measured using annular chambers.

Technical Specification :

Flow rate can be adjusted via valves
Pressure measurement with annular chambers
Investigations on pipes
4 Pipe section of the same lengths and different materials as measured objects
Flow rate measurement using rotameter
Water connections made using quick-release couplings
Pipe sections can be selected via ball valves
Pressure measurement with digital differential pressure manometer, battery-operated.



Laboratory Instrument India

