



Email : sales@laboratoryinstrumentindia.com

Product Name :
Stainless Steel Mesh Sieve Set

Product Code :
LBNY-0004-770011



Description :

Stainless Steel Mesh Sieve Set

Technical Specification :

Stainless Steel Mesh Sieve Set

These economical sieves feature stainless steel mesh and zinc plated steel frames.

The set of nine sieves includes the following mesh sizes:

2.36mm, 1.18mm, 600 micron, 425 micron, 300 micron, 212 micron, 150 micron, 75 micron, and 1.00mm micron

Includes lid and catch pan

Sieve Diameter: 200mm.

Sieve brush – double end brass

Sieve brush – double end soft nylon

Sieve brush – double end hard nylon.



Laboratory Instrument India