




Email : sales@laboratoryinstrumentindia.com

Product Name :
Control Unit For Ventilation System

Product Code :
LBNY-0005-12500052



Description :

Electronic controller with temperature sensor and setting of the reference variable
Quick-release couplings for connection to hot water supply and ventilation system
Temperature control of a ventilation system
3-way mixing valve dn 20
3 thermometers and 1 temperature sensor with display
Circulating pump.

Technical Specification :

Measuring Ranges
Temperature: 1x -50...600°C, 3x 0...120°C
Pump
Max. Flow rate: 60L/min
Power consumption: 70W
Max. head: 4m.



Laboratory Instrument India

