



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Demo Gas Burner Unit

**Product Code :**  
LBNY-0005-12500053



**Description :**

Simulation model of a forced air gas burner and all burner operating states  
Important electrical monitoring points on burner accessible via lab jacks and measuring bridge connectors for fault analysis  
Fault simulation: 9 fault switches in lockable box  
Additional actual system components: flame safeguard, gas pressure and air pressure control devices, ignition transformer, thermostat, safety temperature limiter, gas solenoid valves, fan motor, emergency switch heating system.

**Technical Specification :**

Fully automatic air forced gas burner for natural gas  
Connection load: approx. 160w  
Motor power: 70W  
Power: 15...45kw.



**Laboratory Instrument India**

