



Email : sales@laboratoryinstrumentindia.com

Product Name :
Forced Air Gas Burner Training Panel

Product Code :
LBNY-0005-12500056



Description :

Simulation of the operation of a gas burner
Indication of burner operation using a flickering light
Monitoring of burner operation using flame safeguard
Hazard-free operation of the burner with air
Display of the most important operating pressures on single tube manometers
Diaphragm pump for pressure generation, pumping medium: air.

Technical Specification :

Diaphragm Pump
Max. flow rate: 70L/min at 230mbar
Power consumption: 70W
Gas Burner
Outlet pressure: 0,04bar
Rating class: 13...28kw
Admission pressure: 0,4bar
Measuring Ranges
Pressure: 3x 0...50mbar
Flow rate air: 0,5...5,5m³/h.

Laboratory Instrument India

