



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Oil Tank Safety Trainer

**Product Code :**  
LBNY-0005-12500057



**Description :**

Oil tank according to double-walled  
Oil tank, complete with all necessary filling, outlet and safety devices  
Liquid leak detection system, vacuum leak detection system  
Single-strand oil filter, double-strand oil filter  
Device for overfill protection  
Combination fitting  
Dip stick tube with dip stick and seal  
Tank ventilation  
Pneumatic level gauge  
Filler pipe with seal.

**Technical Specification :**

Oil tank: 350L  
1 leak detector  
1 leak detector with tank.



**Laboratory Instrument India**

