



Email : sales@laboratoryinstrumentindia.com

Product Name :
Oil Tank Safety Trainer

Product Code :
LBNY-0005-12500057



Description :

Oil tank according to double-walled
Oil tank, complete with all necessary filling, outlet and safety devices
Liquid leak detection system, vacuum leak detection system
Single-strand oil filter, double-strand oil filter
Device for overfill protection
Combination fitting
Dip stick tube with dip stick and seal
Tank ventilation
Pneumatic level gauge
Filler pipe with seal.

Technical Specification :

Oil tank: 350L
1 leak detector
1 leak detector with tank.



Laboratory Instrument India

