



Email : sales@laboratoryinstrumentindia.com

Product Name :
LCR Circuit

Product Code :
LBNY-0004-780001



Description :

LCR Circuit, LCR Circuit Equipments, LCR Circuit Tools, LCR Circuit Tool Kits, LCR Circuit Manufacturers, LCR Circuit Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

Technical Specification :

LCR Circuit

This installation is indicative. All equivalent or better solution pedagogically and technically will be accepted :

Compatible with computer interface

Circuit board includes:

Four resistors: 47KOhm, 3.3 KOhm and two 1.0KOhm

Two inductors: 6.8 mH and 2.5 mH

Two capacitors: 3900picoF and 560picoF.

Laboratory Instrument India

