



Email : sales@laboratoryinstrumentindia.com

Product Name :
Instantaneous Gas Heater

Product Code :
LBNY-0005-12500058



Description :

PLC simulation of safety functions
Simulation of 12 faults
Safe operation of heater using air
Simulation of an instantaneous gas heater
Diaphragm pump for pressure generation.

Technical Specification :

Instantaneous heater
Diaphragm pump: max. flow rate: 70L/min
Gas pressure controller: 16...28mbar
Nominal thermal output: 19.2kW
Flow Meters
Water: 75...1000L/h
Air: 0.75...7nm³/h (standard cubic meter per hour).
Manometer
U-tube manometer for air: 10-0-10mbar
Gas: 0...60mbar
Water: 0...6bar.

Laboratory Instrument India

