



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Domestic Gas Supply Training Panel

**Product Code :**  
LBNY-0005-12500059



**Description :**

Exercises on gas pipes

1 gas meter with 1 connection

1 gas meter with 2 connections

2 individually selectable pipe sections

1 pipe section as supply pipe with gas socket and gas connection with ball valve

1 pipe section for finding leaks.

**Technical Specification :**

Gas Pressure Controller

Outlet pressure: 20mbar

Admission pressure: 350mbar

Leak testing device: max. 500mmWC.

Air Pressure Controller

Outlet pressure: 0.01...0.6bar

Max. Inlet pressure: 10bar.



**Laboratory Instrument India**

