




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Three Way Mixing Valve Training Unit

**Product Code :**  
LBNY-0005-12500060



#### **Description :**

Trainer on heating systems and plumbing  
4 bimetallic dial thermometers  
2 water connections for boiler dn 15  
2 heating water connections dn 15  
3-way mixing valve dn 20  
Circulating pump  
1 differential pressure overflow valve  
4 rotameters  
Water connections made using quick-release couplings.

#### **Technical Specification :**

Three-way mixing valve  
Opening time: 100s  
Setting angle: 90°  
Pump  
Max. Flow rate: 60L/min  
Power consumption: 70W  
Max. head: 4m  
Measuring Ranges  
Temperature: 0...100°C  
Flow rate: 3x 150...1700l/h, 1x 40...400l/h  
Dimensions And Weight

---

LxWxH: 1650x700x1850mm  
Weight: approx. 112kg.



**Laboratory Instrument India**