




Email : sales@laboratoryinstrumentindia.com

Product Name :
Three Way Mixing Valve Training Unit

Product Code :
LBNY-0005-12500060



Description :

Trainer on heating systems and plumbing
4 bimetallic dial thermometers
2 water connections for boiler dn 15
2 heating water connections dn 15
3-way mixing valve dn 20
Circulating pump
1 differential pressure overflow valve
4 rotameters
Water connections made using quick-release couplings.

Technical Specification :

Three-way mixing valve
Opening time: 100s
Setting angle: 90°
Pump
Max. Flow rate: 60L/min
Power consumption: 70W
Max. head: 4m
Measuring Ranges
Temperature: 0...100°C
Flow rate: 3x 150...1700l/h, 1x 40...400l/h
Dimensions And Weight

LxWxH: 1650x700x1850mm
Weight: approx. 112kg.



Laboratory Instrument India