



Email : sales@laboratoryinstrumentindia.com

Product Name :
Radiator Training Panel

Product Code :
LBNY-0005-12500061



Description :

Trainer on heating systems and plumbing
2 water connections for boiler DN 15
2 cooling water connections DN 15
4 radiators with thermostatic valve, air bleed and adjustable lock shield valve
Water connections made using quick-release couplings
5 rota meters.

Technical Specification :

Measuring ranges
Flow rate: 1x 0...1000L/h, 4x 0...300L/h
Radiator
Capacity: 3kW
Plate heat exchanger with 10 plates.



Laboratory Instrument India

