



Email : sales@laboratoryinstrumentindia.com

Product Name :
Four Way Mixing Valve Training Panel Unit

Product Code :
LBNY-0005-12500064



Description :

Trainer on heating systems and plumbing
4 rotameters
4 bimetallic dial thermometers
2 water connections for boiler DN 15
4-way mixing valve DN 20
1 differential pressure overflow valve
Circulating pump
2 heating water connections DN 15
Water connections made using quick-release couplings.

Technical Specification :

Pump
Max. Flow rate: 60L/min
Power consumption: 60W
Max. head: 4m
Measuring Ranges
Temperature: 0...100°C
Flow rate: 3x 150...1600l/h, 1x 0...440l/h.
Four-Way Mixing Valve
Opening time: 2min
Setting angle: 90°.



Laboratory Instrument India