



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Temperature Measurement Training Unit

**Product Code :**  
LBNY-0005-12500066



**Description :**

Digital display of the temperatures measured using sensors  
Flow rate measurement with rotameter  
4 different measuring sections, can be selected individually with ball valves  
Water connections made using quick-release couplings  
Investigation of 4 different temperature measuring techniques in the range from 0...60°C.

**Technical Specification :**

Measuring Ranges  
Rota meter: 150...1600L/h  
Thermometer  
Liquid expansion thermometer: 0...60°C  
Resistance thermometer pt100: -50...400°C  
Bimetallic dial thermometer: 0...60°C  
Thermocouple type K: 0...1200°C.



**Laboratory Instrument India**

