



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipe Cleaning Training Panel Unit

Product Code :
LBNY-0005-12500069



Description :

Usage of original plumbing components
Pipe flushing in accordance
Transparent pipe fittings enable optimal observation of the flushing process
Cold water connection using hose
Alternative accessories: flushing compressor (not included)
Contaminants can be introduced into the pipe system via filling openings that can be shut off.

Technical Specification :

Recommended flow pressure: 1...5bar
Automatic mixer, thermostatic mixer
Max. Water temperature: 80C, with temperature
Min. Flow rate: 5L/min

Safety stop
Quick-action bleed valve
Safety valve on the hot water preparation unit
Two handle mixer
Safety fittings
Pipes: Cu / transparent plastic
Simulated hot water preparation
10x hose connection with quick-release coupling
2x cold water connection: Geka coupling 1.



Laboratory Instrument India