



Email : sales@laboratoryinstrumentindia.com

Product Name :
Domestic Heating Circuit Training Panel Unit

Product Code :
LBNY-0005-12500070



Description :

Trainer on heating systems and plumbing
Threaded heater element with switch cabinet
Pressure vessel
1 system of safety devices with pressure reducing valve, shut-off valve, test connections and non-return valve, manometer connection and safety valve
2 spring tube manometers
Thermal discharge safety device
1 bimetallic dial thermometer
2 blow-off lines, 1 filler pipe, 1 drain.

Technical Specification :

Measuring Ranges:
Pressure: 0...10bar
Temperature: 0...120°C
Pressure vessel made of steel
Volume: approx. 5l
Max. operating pressure: 7,8bar
Threaded Heater Element:
Power: 3kw
Protective system: IP 65
Setting range of pressure reducing valve: 1,5...6bar
Thermal discharge safety device: actuated at 95°C.



Laboratory Instrument India