



Email : sales@laboratoryinstrumentindia.com

Product Name :
Protection Of Drinking Water Training Panel Unit

Product Code :
LBNY-0005-12500071



Description :

Devices to prevent the return of contaminated water into the drinking water pipe
Simulation of 3 floors and a cellar of a house
Pipe system made of transparent material
6 outlet valves with hose connections
2 risers with vent
Pipe separator, backflow preventer and riser vent as safety devices
Possible simulations: pipe fracture, pressure drop
3 transparent water tanks
Water meter in feed pipe

Technical Specification :

1 pipe separator: DN 20
2 riser vents: DN 20
2 backflow preventers: DN 20
1 flow rate pipe vent without backflow preventer: DN 20
1 flow rate pipe vent with backflow preventer: DN 20
3 Plexiglass tanks: 6,5L
2 sealing plugs: DN 20
2 pipe sections: DN 20.



Laboratory Instrument India