




Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydronic Balancing Of Radiators Unit

Product Code :
LBNY-0005-12500072



Description :

Expansion vessel, safety group according to overflow valve
5 rotameters
1 surface thermometer for temperature measurement at radiator inlets and outlets
Trainer for heating and plumbing
Water connections using quick-release couplings, diameter inside 13mm, for optional connection to an external heating circuit
6 radiators with thermostatic valve, bleed valve and adjustable lock shield valve.

Technical Specification :

Measuring Ranges
Temperature: 1x -50...300°
Flow rate: 1x 30...320l/h, 4x 50...640l/h
Pump
6 balancing valves: PN 16
Power consumption: 60W
Max. Head: 4m
Max. Flow rate: 60L/min.

Laboratory Instrument India

