



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hydronic Balancing Of Radiators Unit

**Product Code :**  
LBNY-0005-12500072



**Description :**

Expansion vessel, safety group according to overflow valve  
5 rotameters  
1 surface thermometer for temperature measurement at radiator inlets and outlets  
Trainer for heating and plumbing  
Water connections using quick-release couplings, diameter inside 13mm, for optional connection to an external heating circuit  
6 radiators with thermostatic valve, bleed valve and adjustable lock shield valve.

**Technical Specification :**

Measuring Ranges  
Temperature: 1x -50...300°  
Flow rate: 1x 30...320l/h, 4x 50...640l/h  
Pump  
6 balancing valves: PN 16  
Power consumption: 60W  
Max. Head: 4m  
Max. Flow rate: 60L/min.

---

# Laboratory Instrument India

