



Email : sales@laboratoryinstrumentindia.com

Product Name :
Safety Devices Training Unit

Product Code :
LBNY-0005-12500073



Description :

Threaded heater element with switch cabinet
2 spring tube manometers
1 bimetallic dial thermometer
Trainer on heating systems and plumbing
Pressure vessel
2 blow-off lines, 1 filler pipe, 1 drain
Thermal discharge safety device
1 system of safety devices with pressure reducing valve, shut-off valve, test connections and non-return valve, manometer connection and safety valve.

Technical Specification :

Threaded heater element
Power: 3kw
Setting range of pressure reducing valve: 1,5...6bar
Thermal discharge safety device: actuated at 95°C
Pressure vessel made of steel
Max. operating pressure: 7,8bar
Volume: approx. 5l
Measuring ranges
Temperature: 0...120°C
Pressure: 0...10bar.



Laboratory Instrument India