



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Aluminium Trank

**Product Code :**  
LBNY-0004-780003



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Aluminium Trank

**Technical Specification :**

Mechanical bench 01:  
Aluminium trank or equivalent  
L ? 2 m  
Metric scale in mm  
Adjustable end stops  
Level adjustment



**Laboratory Instrument India**