



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sanitation Fittings Training Panel Apparatus

Product Code :
LBNY-0005-12500078



Description :

Investigating the most common sanitation fittings
Possibility to test different fittings simultaneously
Facilities for hot and cold water supply using quick-release couplings
Flushing valve
2 Bourdon tube manometers
2 rotameters
Two handle mixer as water tap
Two handle mixer for a shower
2 adjustable pressure reducing valves
3 dial thermometers.

Technical Specification :

Pressure reducing valve
Outlet pressure: 1,5...6bar, up to 70°C
Max. Initial pressure: 25bar
Flushing valve
DN 20
1,2...5bar
Measuring Ranges
Flow rate: 160...1600L/h
Temperature: 0...100°C
Pressure: 0...6bar



Laboratory Instrument India