



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Sanitation Fittings Training Panel Apparatus

**Product Code :**  
LBNY-0005-12500078



**Description :**

Investigating the most common sanitation fittings  
Possibility to test different fittings simultaneously  
Facilities for hot and cold water supply using quick-release couplings  
Flushing valve  
2 Bourdon tube manometers  
2 rotameters  
Two handle mixer as water tap  
Two handle mixer for a shower  
2 adjustable pressure reducing valves  
3 dial thermometers.

**Technical Specification :**

Pressure reducing valve  
Outlet pressure: 1,5...6bar, up to 70°C  
Max. Initial pressure: 25bar  
Flushing valve  
DN 20  
1,2...5bar  
Measuring Ranges  
Flow rate: 160...1600L/h  
Temperature: 0...100°C  
Pressure: 0...6bar



**Laboratory Instrument India**