



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Refrigerant Filling And Evacuation Equipment

**Product Code :**  
LBNY-0005-12500081



**Description :**

Manometer for intake and high pressure with temperature scale for refrigerant R134a  
High pressure and intake pressure connections  
4-valve assembly aid  
Pulsation-damped manometer for intake pressure, high pressure; vacuum meter  
Portable Equipment  
Vacuum pump and filling balance

**Technical Specification :**

Intake volume: 66L/min, final vacuum: 0,02mbar  
Vacuum Pump: Motor Power: 0,25kw,  
Measuring Ranges  
High Pressure: 0...31bar  
Intake Pressure: -1...8bar  
Vacuum: 0...1000mbar.



**Laboratory Instrument India**

