



Email : sales@laboratoryinstrumentindia.com

Product Name :
Refrigerant Filling And Evacuation Equipment

Product Code :
LBNY-0005-12500081



Description :

Manometer for intake and high pressure with temperature scale for refrigerant R134a
High pressure and intake pressure connections
4-valve assembly aid
Pulsation-damped manometer for intake pressure, high pressure; vacuum meter
Portable Equipment
Vacuum pump and filling balance

Technical Specification :

Intake volume: 66L/min, final vacuum: 0,02mbar
Vacuum Pump: Motor Power: 0,25kw,
Measuring Ranges
High Pressure: 0...31bar
Intake Pressure: -1...8bar
Vacuum: 0...1000mbar.



Laboratory Instrument India

