



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Model Ceiling Air Cooler

Product Code :
LBNY-0005-12500095



Description :

Cutaway models of refrigeration components
Filter/drier, L=127mm, d=52mm, 0,27kg
Ceiling fan air cooler, 440x430x130mm, 8,6kg
Thermostatic expansion valve, 79x32x65mm

Technical Specification :

Expansion vessel with internal pressure compensation
Ceiling fan air cooler
Capacity: 0,33kw
Cooling surface area: 2,2m²
Volumetric air flow rate: 260m³/h.



Laboratory Instrument India