



Email : sales@laboratoryinstrumentindia.com

Product Name :
Replacement Of Refrigeration Components Apparatus

Product Code :
LBNY-0005-12500097



Description :

Complete refrigeration circuit with open compressor, air-cooled condenser, receiver and air cooling evaporator
Replaceable components: compressor, filter/drier, pressure switch and solenoid valve
Trainer from the practical series for the training of mechatronics engineers for refrigeration
Recommended accessories evacuation equipment.
Service exercises in refrigeration systems
Overall arrangement on robust workbench.

Technical Specification :

Open refrigerant compressor
Maximum system pressure: 18bar
Refrigeration capacity: 1245W at -5°C/36°C
Drive power: 0,75kw
Pressure switch
HP: 10...32bar
LP: -0.2...7bar.



Laboratory Instrument India

