



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Replacement Of Refrigeration Components Apparatus

**Product Code :**  
LBNY-0005-12500097



**Description :**

Complete refrigeration circuit with open compressor, air-cooled condenser, receiver and air cooling evaporator  
Replaceable components: compressor, filter/drier, pressure switch and solenoid valve  
Trainer from the practical series for the training of mechatronics engineers for refrigeration  
Recommended accessories evacuation equipment.  
Service exercises in refrigeration systems  
Overall arrangement on robust workbench.

**Technical Specification :**

Open refrigerant compressor  
Maximum system pressure: 18bar  
Refrigeration capacity: 1245W at -5°C/36°C  
Drive power: 0,75kw  
Pressure switch  
HP: 10...32bar  
LP: -0.2...7bar.



**Laboratory Instrument India**

