



Email : sales@laboratoryinstrumentindia.com

Product Name :
Refrigerator Unit Model

Product Code :
LBNY-0005-125000101



Description :

Model of a refrigerator to plug onto the base unit
Electric heater to generate the cooling load
Expansion elements selectable via solenoid valves: expansion valve or capillary tube
Training system with hsi technology
Chamber with transparent front
Sensors to record temperature and pressure
Refrigeration chamber with evaporator, fan and cooling load
Operation of solenoid valves, fan, heater and fault simulation via software.

Technical Specification :

Refrigeration chamber, LxWxH: 270x270x220mm
Capillary tube: length 2m
Electric PTC heater as cooling load: 210W
Measuring Ranges:
Pressure: -1...9bar
Temperature: 3x -50...50°C.

Laboratory Instrument India

