



Email : sales@laboratoryinstrumentindia.com

Product Name :
Refrigeration Components For Experimentsig

Product Code :
LBNY-0005-125000102



Description :

Components for the setup of simple refrigeration circuits in conjunction
Pressure-controlled expansion valve
Thermostatic expansion valve
2 thermostats with different temperature ranges
Heat exchanger as super heater
3-pin circuit breaker
Sight glass with humidity indicator
Components mounted on plates ready to connect
Abrasion-proof symbols and labels on plates
Assembly aid: 4-port directional valve with sight glass.

Technical Specification :

Thermostatic expansion valve
Valve insert rated capacity: 0,3kW
Evaporation temperature: -45°...20°c
Circuit breaker, 3-pin: 3kW
Pressure-controlled expansion valve
Pressure range: 1...7bar
2 electric thermostats
25...5°C
5...25°C
Heat Exchanger

Capacity: approx. 2.5kW at -10°C evaporation temperature.



Laboratory Instrument India