



Email : sales@laboratoryinstrumentindia.com

Product Name :
Refrigeration System With Open Compressor

Product Code :
LBNY-0005-125000106



Description :

Open Compressor
Condenser with fan transfer area: 2,5m
Capacity: approx. 1935W at 25C air temperature (ambient) / t=15C
Refrigeration capacity: approx. 845w (at speed: 975min-1 and -10/40c)
Speed: 500...1000min-1

Technical Specification :

Evaporator
Heater power: 500W
Transfer area: 3,62m
Capacity: 460w at 3c air temperature (chamber) / t=13c
Measuring Ranges:
Flow rate r134a: 1x 2...15g/s
Power: 1x 0...500W
Pressure: 2x -1...15bar
Speed: 1x 0...1000min-1
Temperature: 4x -5...105c, 1x -50...250c.

Laboratory Instrument India

