



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Refrigeration System With Open Compressor

**Product Code :**  
LBNY-0005-125000106



**Description :**

Open Compressor  
Condenser with fan transfer area: 2,5m  
Capacity: approx. 1935W at 25C air temperature (ambient) / t=15C  
Refrigeration capacity: approx. 845w (at speed: 975min-1 and -10/40c)  
Speed: 500...1000min-1

**Technical Specification :**

Evaporator  
Heater power: 500W  
Transfer area: 3,62m  
Capacity: 460w at 3c air temperature (chamber) / t=13c  
Measuring Ranges:  
Flow rate r134a: 1x 2...15g/s  
Power: 1x 0...500W  
Pressure: 2x -1...15bar  
Speed: 1x 0...1000min-1  
Temperature: 4x -5...105c, 1x -50...250c.

---

# Laboratory Instrument India

