



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pulley With Walnut of Fixing On The Bench

Product Code :
LBNY-0004-780006



Description :

Pulley With Walnut of Fixing On The Bench

Technical Specification :

Mechanical bench accessory 04:

Pulley with walnut of fixing on the bench

Pulley Bench Mounting which consists of a metal pulley 70mm / 2 1/2" dia, mounted in a metal frame which can be clamped on to a table upto 42mm / 1 5/8" thick.

The frame also has a hole for a rod of upto 12mm diameter to pass through, if wanted the pulley can be mounted on a rod.



Laboratory Instrument India