



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Pulley With Walnut of Fixing On The Bench

**Product Code :**  
LBNY-0004-780006



**Description :**

Pulley With Walnut of Fixing On The Bench

**Technical Specification :**

Mechanical bench accessory 04:

Pulley with walnut of fixing on the bench

Pulley Bench Mounting which consists of a metal pulley 70mm / 2 1/2" dia, mounted in a metal frame which can be clamped on to a table upto 42mm / 1 5/8" thick.

The frame also has a hole for a rod of upto 12mm diameter to pass through, if wanted the pulley can be mounted on a rod.



**Laboratory Instrument India**