




Email : sales@laboratoryinstrumentindia.com

Product Name :
Water Cooled Four Stroke Diesel Engine

Product Code :
LBNY-0005-125000118



Description :

Engine complete with fuel hose and exhaust gas temperature sensor
Cooling water circuit with circulating pump, flow meter and temperature sensors
Water-cooled single-cylinder four-stroke diesel engine with swirl chamber for installation .
Force transmission to brake via elastic claw coupling
Fuel hose with self-sealing quick-release coupling
Engine mounted on base plate.

Technical Specification :

Single-cylinder diesel engine with indirect injection
Power output: 5,1kW at 3000min-1
Bore: 75mm, stroke: 70mm
Cooling water circuit
Pump: max. 640L/h, max. 100m
Capacity: 2,5L
Measuring ranges
Flow rate: 30...300L/h
Temperature: 2x 0...100°C.

Laboratory Instrument India

