



Email : sales@laboratoryinstrumentindia.com

Product Name :
Air Compressor Test Unit

Product Code :
LBNY-0005-125000121



Description :

Pressure vessel with valve for draining
Display at switch cabinet to indicate active power
Test unit for a small industrial compressor
1-cylinder compressor

Technical Specification :

Drive motor
Speed: 1410min-1
Power consumption: 0,25kw
Pressure vessel
Max. Pressure load: 10bar
Capacity: 10L
Compressor
Intake capacity: 115L/min
Power consumption: 0,7kw
Max. Pressure: 10bar
Measuring ranges
Stopwatch: 0...9h59min59s 99/100
Manometer: 0...16bar.



Laboratory Instrument India