



Email : sales@laboratoryinstrumentindia.com

Product Name :
Simulation Of A Two Stage Air Compressor

Product Code :
LBNY-0005-125000124



Description :

9 potentiometers for the setting of
Relative humidity: 0...100%
Coolant flow rate: 0...100kg/h at 15°C
Intake pressure: 0...2bar abs.
Intake temperature: 0...100°C
Flow control valve setting: 0...100%
Vessel, nominal pressure: 0...50bar
Dead space: 0...100%
Motor speed: 0...1000min⁻¹
Vessel capacity: 0...1000l
Inputs and outputs
4x digital in/out
16x analogue in, 1x analogue out.

Technical Specification :

Simulation Of A Two Stage Air Compressor

Laboratory Instrument India

