



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Parallel Flow Module Apparatus

**Product Code :**  
LBNY-0005-125000126



**Description :**

2 thermocouples type K: Temperature measurement at the top of the rod-type heater and at the inner surface of the ring heater.  
Mounted with quick-action fasteners.

**Technical Specification :**

Air duct: D=60mm  
Ring heater  
Length: 30mm  
Power: 220W  
Diameter: 60mm  
Rod-Type Heater  
Length: 130mm  
Power: 250W  
Diameter: 8mm



**Laboratory Instrument India**

