



Email : sales@laboratoryinstrumentindia.com

Product Name :
Parallel Flow Module Apparatus

Product Code :
LBNY-0005-125000126



Description :

2 thermocouples type K: Temperature measurement at the top of the rod-type heater and at the inner surface of the ring heater.
Mounted with quick-action fasteners.

Technical Specification :

Air duct: D=60mm
Ring heater
Length: 30mm
Power: 220W
Diameter: 60mm
Rod-Type Heater
Length: 130mm
Power: 250W
Diameter: 8mm



Laboratory Instrument India

