



Email : sales@laboratoryinstrumentindia.com

Product Name :
Set of 12 Harmonic Springs With Different Stiffness

Product Code :
LBNY-0004-780008



Description :

Set of 12 Harmonic Springs With Different Stiffness

Technical Specification :

Mechanical bench accessory 06:

Set of 12 harmonic springs with different stiffness.

12 flexible extension springs, made out of quality steel

Looped ends make this extension spring an excellent tool for studying motion, mass and potential energy.

Perfect for physics experimentation in the classroom or at home.

Sturdy construction allows for long time use and handling.



Laboratory Instrument India