



Email : sales@laboratoryinstrumentindia.com

Product Name :
Heat Transfer Bench Apparatus

Product Code :
LBNY-0005-125000134



Description :

Display and control unit indicates air temperature, heater temperature and heating power
2 heaters can be inserted in a bundle or alone
Air flow rate adjustable
Convective heat transfer
Exchangeable pipe bundle as heat exchanger
Pitot tube and pressure measurement station to determine air velocity.

Technical Specification :

Duct
Length: 1540mm
Cross-section: 150x150mm
Pipe bundle
23 pipes (D=13mm)
23 pipes (D=10mm)
Fan
Max. flow rate: 2160m³/h
Power: 1,5kw
2 heaters
Power: 1x250W (D=13mm)
Power: 1x220w (D=10mm)



Laboratory Instrument India