




Email : sales@laboratoryinstrumentindia.com

Product Name :
MOT Control of a DC Motor

Product Code :
LBNY-0005-1260003



Description :

Industrial type components must be educationally adapted by using a modular panel system to permit the step by step assembling from the simplest circuit to the most complex system. This trainer has a modular structure and it consists of panels installed on a vertical frame. It is supplied with a theoretical and practical manual. This laboratory is designed for the study of Control of a DC motor to allow the student a practical training, based on the performance of guided experiments. The modularity of this system can give the students a direct and immediate approach to the topics, offering the opportunity to study various subjects, performing several experiments as following:

P, PI and PID controls of the speed of a DC motor using the CHR method
Control of a DC generator.

Technical Specification :

MOT Control of a DC Motor



Laboratory Instrument India