




Email : sales@laboratoryinstrumentindia.com

Product Name :
Flow Control

Product Code :
LBNY-0005-1260004



Description :

This trainer has a modular structure and it consists of panels installed on a vertical frame. It is supplied with a theoretical and practical manual. Industrial type components are educationally adapted by using a modular panel system to permit the step by step assembling from the simplest circuit to the most complex system. This laboratory is designed for the study of Flow control application to allow the student a practical training, based on the performance of guided experiments. The modularity of this system can give the students a direct and immediate approach to the topics, offering the opportunity to study various subjects, performing several experiments as following:

Analysis of flow control features
Measurement systems features.

Technical Specification :

Flow Control



Laboratory Instrument India