



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Bench with Transparent Components

Product Code :
LBNY-0005-1260005



Description :

This Trainer Is Composed As Follows:

Work bench complete with:

4 wheels (the 2 front wheels are lockable)

3 drawers with locking systems with key for storing components.

50mm (horizontally and vertically) for quick clamping of hydraulic components without using additional tools.

1 clipboard located on the side of the work bench.

Metal working panel fully plastic coated; panel dimensions:950 x 1050 mm. With 5 mm holes every.

Technical Specification :

Hydraulic Bench with Transparent Components



Laboratory Instrument India