




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Continuous Automatic Control

**Product Code :**  
LBNY-0005-1260007



**Description :**

This trainer has a modular structure and it consists of panels installed on a vertical frame. The analysis of the interactions between controllers and processes is complicated by the possible presence of noise sometimes, the latter can trigger a series of oscillations with consequences, potentially serious, for the process. This laboratory is designed to study processes and controllers suitably combined to simulate and to study the most common problems related to the Continuous Automatic Control. In this configuration, the student can study the causes of the above mentioned instabilities, in order to find possible solutions.

**Technical Specification :**

Continuous Automatic Control



**Laboratory Instrument India**