



Email : sales@laboratoryinstrumentindia.com

Product Name :
Single Phase Electro Compressor

Product Code :
LBNY-0005-1260008



Description :

With very low noise level for small capacity.
Air intake: 50 litres/min.
Max. working pressure: 8 bar or 116 PSI
Tank capacity: 24 litres
Noise level: 40 dB
Motor power: 0.34 kW or 0.45 HP
Dimensions: 40 x 40 x 60 (h) cm.
Weight: 25 kg.

Technical Specification :

Single Phase Electro Compressor



Laboratory Instrument India