




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Level Control

**Product Code :**  
LBNY-0005-1260009



#### **Description :**

This trainer has a modular structure and it consists of panels installed on a vertical frame. It is supplied with a theoretical and practical manual. Industrial type components are educationally adapted by using a modular panel system to permit the step by step assembling from the simplest circuit to the most complex system. This laboratory is designed for the study of the Level Control Application to allow the student a practical training, based on the performance of guided experiments. The modularity of this system can give the students a direct and immediate approach to the topics, offering the opportunity to study various subjects, performing several experiments as following:

#### **Pumps features:**

Measurement systems features  
Level control with adjuster P  
Level control with adjuster PI  
Level process features  
Level control with adjuster PID  
Level control with a 2 position controller.

#### **Technical Specification :**

Level Control

---

# Laboratory Instrument India

