




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Light Control

**Product Code :**  
LBNY-0005-12600010



**Description :**

This trainer has a modular structure and it consists of panels installed on a vertical frame. It is supplied with a theoretical and practical manual. This laboratory is designed for the study of the Light Control application to allow the student a practical training, based on the performance of guided experiments. Industrial type components are educationally adapted by using a modular panel system to permit the step by step assembling from the simplest circuit to the most complex system. The modularity of this system can give the students a direct and immediate approach to the topics, offering the opportunity to study various subjects, performing several experiments as following:

P, PI and PID light controls using the CHR method.

**Technical Specification :**

Light Control



**Laboratory Instrument India**