




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Temperature Control

**Product Code :**  
LBNY-0005-12600012



**Description :**

Industrial type components are educationally adapted by using a modular panel system to permit the step by step assembling from the simplest circuit to the most complex system. This laboratory is designed for the study of the temperature control application to allow the student a practical training, based on the performance of guided experiments.

With this laboratory it is possible to perform the following experiments:

The two position controller with delayed feedback in the temperature process

The two position controller with elastic feedback in the temperature process

The two position controller in the temperature process

P, PI and PID controls of the temperature process using the CHR method

The three range controller in the temperature process.

**Technical Specification :**

Temperature Control



**Laboratory Instrument India**

